



Procalcitonin (PCT)

Safely reduce antibiotic exposure

B·R·A·H·M·S PCT:

An effective tool for antibiotic stewardship



The challenge

Antibiotic resistance – an increasing threat to public health


Antibiotics (ABx) are a limited resource. At the current pace of injudicious use, all antibiotics will soon become ineffective. The **WHO Global Action Plan on Antimicrobial Resistance**, 2015, emphasizes that antimicrobial resistance is a crisis that must be managed with the utmost urgency.¹

The **Interagency Coordination Group (ICG) on Antimicrobial Resistance** reported in 2019 that, “unless the world acts urgently, antimicrobial resistance will have disastrous impact within a generation.” Deaths due to drug-resistant diseases “could increase to 10 million deaths globally per year by 2050.”²

1/3 antibiotic prescriptions are unnecessary³



How does resistance to antibiotics develop?



A potential for change

Thermo Scientific™ B·R·A·H·M·S PCT™ supports responsible use of antibiotics to prolong their effectiveness


“For adults with an initial diagnosis of sepsis or septic shock and adequate source control where optimal duration of therapy is unclear, we suggest **using procalcitonin AND clinical evaluation to decide when to discontinue antimicrobials** over clinical evaluation alone.”

Surviving Sepsis Campaign
Guideline 2021⁴

Procalcitonin is the **ONLY** recommended biomarker for antibiotic stewardship in sepsis and LRTI.

World Health Organization
Essential In Vitro Diagnostics List 2019⁵


B·R·A·H·M·S PCT-aided antibiotic therapy has the potential to ...



Use of B·R·A·H·M·S PCT reduces antibiotic exposure

Strong evidence supports safe reduction of antibiotics using PCT-aided antibiotic stewardship protocols

- Proven utility across diverse clinical settings
- Reproducible, randomized clinical trials with evidence spanning over 20 years




Proven efficacy:
-16% to -74%
 antibiotic exposure

No
 adverse impact
 on outcome

Figure 1. Relative reduction in ABx exposure with PCT-aided therapy
 ABx exposure in control group is normalized to 100, shown by a blue bar at the top. The gray bar depicts the relative exposure in PCT group and the orange bar shows the relative ABx exposure reduction. All studies reported significant reduction in ABx exposure.

- ABx exposure in control group (normalized to 100)
- Relative ABx exposure in PCT-aided group
- Relative ABx reduction


Benefits of reduced ABx exposure achieved by PCT-aided antibiotic therapy



PCT-aided antibiotic therapy led to a reduction in overall treatment costs: 26% in sepsis patients and 18% in LRTI patients²¹

Safe approach


B·R·A·H·M·S PCT-aided reduction in antibiotic exposure has no negative impact on survival



Overall better survival in the PCT group
Odds ratio: 0.89


Figure 2. Effect of PCT-guided antibiotic treatment on clinical outcomes in ICU patients – results of a patient level meta-analysis of 11 randomized trials (2252 PCT group patients, 2230 control group patients)²²

Lower 28-day mortality




PCT-aided ABx therapy contributes to improved clinical outcomes

Explanation for the 28-day survival benefit seen in the PROGRESS study¹²



* Reduction of diarrhea leads to lower incidence of severe and life-threatening consequences like electrolyte disturbances, dehydration, cardiovascular instability, and acute kidney injury.
MDRO Multidrug-resistant organisms

When to start antibiotics?



PCT values should always be interpreted in context of the patient's clinical condition.

PCT levels increase 3–6 hours after bacterial challenge and return to normal as the infection is resolved (Figure 3)^{24,25,26}

- High specificity and sensitivity for bacterial infection
- Indicator for disease severity and treatment response





Figure 3. Kinetics of PCT^{24,25}

When to stop antibiotics?



PCT values should always be interpreted in context of the patient's clinical condition.

Daily monitoring of PCT course allows for **customized ABx treatment duration, hence reduced ABx exposure**



Ensure using the quality assay for **SAFE** clinical decision making


PCT cut-offs and clinical algorithms were established by use of the global reference standard Thermo Scientific B-R-A-H-M-S PCT sensitive KRYPTOR™ assay and are valid solely for B-R-A-H-M-S PCT assays.

Efficacy of PCT in adults with LRTI symptoms

Patients in the ED

Is it bacterial infection?

As much as 75% of all antibiotic doses are prescribed for acute respiratory tract infections, despite their mainly viral cause.⁷ PCT-aided therapy in such patients allows reduction in ABx exposure without any adverse impact on outcome.⁶



Data from: Effect of Procalcitonin-Based Guidelines vs Standard Guidelines on Antibiotic Use in Lower Respiratory Tract Infections (ProHOSP)⁶

Largest prospective, multicenter, randomized controlled trial with PCT in LRTI patients presenting to EDs:

- 1359 LRTI patients, 6 centers
- PCT group (n=671), control group (n=688)

* % reduction related to non PCT-aided group



0.25 ng/ml

B·R·A·H·M·S PCT
LRTI cut-off

Community-acquired pneumonia (CAP)

Tailor the treatment duration in hospitalized patients





Figure 4. Meta-analysis data for 2027 patients hospitalized for CAP, total exposure of ABx in median days: PCT group = 6 days, control group = 10 days²⁷

Acute COPD exacerbations

Does every exacerbation require ABx?

- Significant sustained reduction in total antibiotic exposure for up to 6 months⁹
- No decrease in mean time to next exacerbation⁹
- No increase in lung function decline⁹



If it is viral, antibiotics will not help. PCT testing can aid in decision making on antibiotic therapy.

Efficacy of PCT in adults in Intensive Care Units

How to know the appropriateness of an empiric antibiotic?

Effective antibiotic treatment is reflected by declining PCT values,²⁸ consistent with its half-life time of about 20–24 hours.²⁵ Serial determinations of PCT can be used to monitor the course of infection in sepsis patients. Appropriate empiric antibiotic therapy is associated with a significant decline in PCT.²⁸





Figure 5. Typical course of PCT serum level according to patient's response to antibiotic treatment (n=109)²⁹

0.5 ng/ml

B·R·A·H·M·S PCT
sepsis cut-off



Efficacy and safety in critically ill patients




Figure 6. Median duration of ABx treatment in PCT-guided group = 5 days, in control group = 7 days¹³

Data from: The Stop Antibiotics on Procalcitonin Guidance Study¹³

- Largest prospective, multicenter, randomized, controlled, open-label intervention trial with PCT in critically ill patients
- Conducted in the Netherlands – a healthcare system with comparatively low use of ABx²³
- 1575 critically ill patients, 15 centers

Multi-disciplinary approach to PCT-aided ABS

Several recent studies have emphasized the importance of a multi-disciplinary approach for successfully leveraging the benefit of PCT.^{30,31,32} A pharmacist led PCT-protocol implementation study showed not only a reduction in days of ABx therapy from 9 days to 6 days but also improved significantly the appropriate ordering of PCT tests, subsequently resulting in higher antibiotic discontinuation rates (Figure 7).³²





Figure 7. Impact of pharmacist-led PCT protocol on appropriateness of ordering of PCT levels and its impact on antibiotic prescriptions³²

B·R·A·H·M·S PCT

PCT-aided antibiotic stewardship protocols lead to



The economic impact of PCT-guided treatment has been studied through health economic modeling in various settings:

- Sepsis patients – ICU³³
- Acute Respiratory Infections – inpatient, ICU, outpatient³⁴
- COPD exacerbation – inpatient³⁵

Treatment cost reductions ranging from 9% to 49% have been demonstrated across various countries.^{21,33,36,37}

The cost of testing for PCT is more than offset by downstream cost savings.




“PCT helps me to prescribe antibiotics rationally and thus to save their power for future generations.”

References

1. WHO Global Action Plan Antimicrobial Resistance, 2015
2. Interagency Coordination Group (IACG) on Antimicrobial Resistance. Report to the Secretary-General of the United Nations, April 2019. https://www.who.int/antimicrobial-resistance/interagency-coordination-group/IACG_final_summary_EN.pdf. Date of access: 10.03.2022
3. Centers for Disease Control and Prevention (CDC), Press Release, May 3, 2016: <http://www.cdc.gov/media/releases/2016/p0503-unnecessary-prescriptions.html>. Date of access: 10.03.2022
4. Evans et al., *Int Care Med* 2021; 47(11): 1181-1247. doi: 10.1007/s00134-021-06506-y
5. World Health Organization (2021). The selection and use of essential in vitro diagnostics: report of the third meeting of the WHO Strategic Advisory Group of Experts on In Vitro Diagnostics, 2020 (including the third WHO model list of essential in vitro diagnostics). Geneva: World Health Organization; 2021 (WHO Technical Report Series, No. 1031). License: CC BY NC-SA 3.0 IGO. Available at: <https://www.who.int/publications/i/item/9789240019102/>. Date of access: 10.03.2022
6. Schuetz et al., *JAMA* 2009; 302 (10): 1059-1066. doi: 10.1001/jama.2009.1297
7. Christ-Crain et al., *Lancet* 2004; 363 (9409): 600-607. doi: 10.1016/S0140-6736(04)15591-8
8. Christ-Crain et al., *Am J Respir Crit Care Med* 2006; 174 (1): 84-93. doi: 10.1164/rccm.200512-19220C
9. Stolz et al., *Chest* 2007; 131 (1): 9-19. doi: 10.1378/chest.06-1500
10. Long et al., *Respirology* 2011; 16 (5): 819-824. doi: 10.1111/j.1440-1843.2011.01978.x
11. Briel et al., *Arch Intern Med* 2008; 168 (18): 2000-2007. doi: 10.1001/archinte.168.18.2000
12. Kyriazopoulou et al., *Am J Respir Crit Care Med* 2021; 203 (2): 202-210. doi: 10.1164/rccm.202004-12010C
13. de Jong et al., *Lancet Infect Dis* 2016; 3099: 1-9
14. Bouadma et al., *Lancet* 2010; 375 (9713): 463-474. doi: 10.1016/S0140-6736(09)61879-1
15. Nobre et al., *Am J Respir Crit Care Med* 2008; 177 (5): 498-505. doi: 10.1164/rccm.200708-12380C
16. Stolz et al., *Eur Respir J* 2009; 34: 1364-1375. doi: 10.1183/09031936.00053209
17. Hochreiter et al., *Crit Care* 2009; 13 (3): R83. doi: 10.1186/cc7903
18. Schroeder et al., *Langenbecks Arch Surg* 2009; 394: 221-226. doi: 10.1007/s00423-008-0432-1
19. Schuetz et al., *The Lancet* 2018, 18 (1): P95-107. doi: 10.1016/S1473-3099(17)30592-3
20. Broyles, *Open Forum Infect Dis* 2017; 4 (4): ofx213. doi: 10.1093/ofid/ofx213
21. Mewes JC et al., *PLoS One* 2019; 14 (4): e0214222. doi: 10.1371/journal.pone.0214222
22. Wirz et al., *Crit Care* 2018; 22: 191. doi: 10.1186/s13054-018-2125-7
23. www.ecdc.europa.eu/en/antimicrobial-consumption/database/rates-country. Date of access: 10.03.2022
24. Harbarth et al., *Am J Respir Crit Care Med* 2001; 164 (3): 396-402. doi: 10.1164/ajrccm.164.3.2009052
25. Meisner M, *Procalcitonin – Biochemistry and Clinical Diagnosis*, ISBN 978-3-8374-1241-3, UNI-MED, Bremen 2010
26. Müller et al., *Crit Care Med* 2000, 28 (4): 977-983. doi: 10.1097/00003246-200004000-00011
27. Schuetz et al., *Clin Infect Dis* 2012; 55 (5): 651-62. doi: 10.1093/cid/cis464
28. Charles et al., *Crit Care* 2009; 13 (2): R38. doi: 10.1186/cc7751
29. Stüber, 21st Int. Congress of Intensive Care and Emergency Medicine (ISICEM), Brussels 2001
30. Willmon et al., *Hospital Pharmacy* 2021; 56 (5): 501-506. doi: 10.1177/0018578720925385
31. Teehan et al., *Hospital Pharmacy* 2021; 56 (5): 569-575. doi: 10.1177/0018578720931450
32. Cowey & Wilson, *Am J Health Syst Pharm* 2022; 79 (4): 239-243. doi: 10.1093/ajhp/zxab413
33. Kip et al., *J Med Econ* 2015; 4: 1-10. doi: 10.3389/fpubh.2016.00115
34. Schuetz et al., *Clin Chem Lab Med* 2015; 53(4): 583-592. doi: 10.1515/cclm-2014-1015
35. Maas et al., Poster, ISPOR 18th Annual European Congress, Milan, 2015
36. Steuten, Cost effectiveness of biomarker testing in sepsis patients; ISICEM Brussels, March 15-19, 2015
37. Voermans et al., *OMICS A Journal of Integrative Biology* 2019; 23(10): 508-515. doi: 10.1089%2Fomi.2019.0113

B·R·A·H·M·S PCT

Secured clinical decision making independent of platform



B·R·A·H·M·S PCT sensitive KRYPTOR



ALINITY i B·R·A·H·M·S PCT
ARCHITECT B·R·A·H·M·S PCT



VIDAS® B·R·A·H·M·S PCT



LIAISON® B·R·A·H·M·S PCT II GEN



Lumipulse® G B·R·A·H·M·S PCT



PATHFAST™ B·R·A·H·M·S PCT



VITROS® B·R·A·H·M·S PCT



ELECSYS® B·R·A·H·M·S PCT



ADVIA Centaur® B·R·A·H·M·S PCT
Atellica® IM B·R·A·H·M·S PCT
Dimension® EXL™ LOCI® B·R·A·H·M·S PCT

Not all intended uses and applications mentioned in this printing are registered for every product in every country.

Clinical Diagnostics

Germany

Thermo Fisher Scientific
B·R·A·H·M·S GmbH
Neuendorfstr. 25
16761 Hennigsdorf,
Germany

+49 (0)3302 883 0
+49 (0)3302 883 100 fax
info.pct@thermofisher.com
www.thermoscientific.com/brahms

USA

Thermo Fisher Scientific
8365 Valley Pike
Middletown, VA 22645

(800) 232 3342
customerservice.diagnostics.mtn@thermofisher.com
www.thermofisher.com/procalcitonin

Learn more at thermoscientific.com/procalcitonin
or email us at info.pct@thermofisher.com



Not all products are CE marked or have 510(k) clearance for sale in the U.S. Availability of products in each country depends on local regulatory marketing authorization status.

© 2023 Thermo Fisher Scientific Inc. All rights reserved. B·R·A·H·M·S PCT and all other trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. ADVIA Centaur, Atellica IM, Dimension EXL and LOCI are registered and protected trademarks belonging to Siemens Healthcare Diagnostics. ADVIA Centaur B·R·A·H·M·S PCT, Atellica IM B·R·A·H·M·S PCT, and Dimension EXL LOCI B·R·A·H·M·S PCT are products of Siemens Healthcare Diagnostics licensed from Thermo Fisher Scientific. ALINITY i and ARCHITECT B·R·A·H·M·S PCT are products of Abbott licensed from Thermo Fisher Scientific. Elecsys is a registered and protected trademark belonging to Roche or one of its subsidiaries. Elecsys B·R·A·H·M·S PCT is a product of Roche licensed from Thermo Fisher Scientific. LIAISON is a registered and protected trademark belonging to DiaSorin S.p.A. LIAISON B·R·A·H·M·S PCT II GEN is a product of DiaSorin S.p.A. licensed from Thermo Fisher Scientific. Lumipulse is a registered trademark of Fujirebio Inc. in Japan and in other countries. Lumipulse G B·R·A·H·M·S PCT is a product of Fujirebio Inc. licensed from Thermo Fisher Scientific. PATHFAST is a registered trademark of LSI Medience Corporation and manufactured in Japan. PATHFAST (reagents, accessories and the device) is a product of PHC Europe B.V. and sold via distributors in EMEA. PATHFAST B·R·A·H·M·S PCT is a product of LSI Medience Corporation licensed from Thermo Fisher Scientific. VIDAS is a registered trademark of bioMérieux S.A. or one of its subsidiaries. VIDAS B·R·A·H·M·S PCT is a product of bioMérieux licensed from Thermo Fisher Scientific. VITROS is a trademark of Ortho Clinical Diagnostics. VITROS B·R·A·H·M·S PCT is a product of Ortho Clinical Diagnostics licensed from Thermo Fisher Scientific. Ortho Clinical Diagnostics is a subsidiary of QuidelOrtho Corporation. KRYPTOR is a trademark of Cisbio Bioassays, licensed for use by B·R·A·H·M·S GmbH, a part of Thermo Fisher Scientific. Patents: www.brahms.de/patents. **BMKT000612.1**